

Amplification of Nucleic Acids

Jennifer L. Hunt, MD, MEd

Aubrey J. Hough Jr, MD, Endowed Professor of Pathology

Chair of Pathology and Laboratory Medicine

University of Arkansas for Medical Sciences

jhunt2@uams.edu

Analysis of DNA

Polymerase Chain Reaction

Polymerase Chain Reaction

- Discovered in 1983 (Kary Mullis, Cetus)
- Patented technique
- Patent sold to Roche in 1991 for \$300 Million




Polymerase Chain Reaction

- Problems before PCR
 - Quantity of DNA was small
 - The target was a tiny fraction of all DNA



Quiz



Polymerase Chain Reaction


- PCR enabled
 - Amplification for quantity
 - Target specific enrichment


Polymerase Chain Reaction

- Template DNA
- Primers
- Taq enzyme
- Nucleotide building blocks






Template DNA: Bases





Double Stranded DNA




Double Stranded DNA





Primer Design

- Primers
 - Short DNA sequence (18-24 nucleotides)
 - Forward and Reverse
 - “Sense & Antisense”

Primer Design



Primer Binding




Polymerase

- Taq Enzyme
 - *Thermus Aquaticus*
 - DNA polymerase
 - Heat activated

Nucleotides

	Base name	Nucleotide name	
A	Adenine	Deoxyadenosine	Purines
G	Guanine	Deoxyguanosine	
C	Cytosine	Deoxycytidine	Pyrimidines
T	Thymine	Deoxythymidine	



PCR Amplification

- Occurs in a machine: Thermocycler (thermal cycler)
- Series of 3 different temperatures = 1 Cycle




Cycle 1

Cycle 2

Denaturation Annealing Extension

Denaturation Annealing




90

72

55



PCR Amplification




PCR Product Analysis

- What do we do after PCR?
 - Detect the sequence in the product
 - The actual sequence (*Sequencing*)
 - Differences that cause altered migration (*MSI*)
 - Detect the amount of product
 - Quantitative (*Q-PCR*)
 - Semi-quantitative (*LOH*)


PCR product detection

- Gel: Agarose or polyacrylamide
- Capillary electrophoresis
- Quantitative PCR


Agarose Gel Electrophoresis





Capillary Electrophoresis




Capillary Electrophoresis




Capillary Gel Electrophoresis






Size of PCR product (basepairs)






Relative amount of PCR product




RNA: Bases (A, C, G, U)




Gene Structure and Translation


DNA




mRNA



Protein



Reverse Transcription – PCR



Reverse Transcription - PCR



cDNA

RT-PCR Product Analysis

- What do we do after RT-PCR?
 - Detect the presence or absence of product
 - Qualitative
 - Detect the sequence in the product, which represents the sequence of the DNA
 - Translocations
 - Detect the amount of product, which represents the amount of expression of a gene
 - Quantitative (*Q-RT-PCR*)



GGCTTT CGG AGATGTTT GATAG CGACGGG AATTTTAACTTTCTCACCTTCTGGGATC

GGCTTT CGG AGATGTTT GATAG CGACGGG AATTTTAACTTTCTCACCTTCTGGGATC

GGCTTT CGG AGATGTTT GATAG CGACGGG AATTTTAACTTTCTCACCTTCTGGGATC